

# LH-133 Series ▲ Customize your limb holder

Check the desired characteristics to create your own limb holder.

Transcribe the corresponding letters and numbers to create your product code number.

LH \_\_\_\_\_ -133 \_\_\_\_\_  
          **A**          **B**          **C**          **D**          **E**

## A ▶ Cuff style

\_\_\_  Regular; if you choose this style, **leave blank** and go to **B**

C  Contoured

CF  Contoured with felt liner



Regular



Contoured

## B ▶ Strap

1  One strap

\_\_\_  Two straps; if you choose this attachment, **leave blank** and go to **C**

## C ▶ Cuff attachment

\_\_\_  Standard; if you choose this cuff, **leave blank** and go to **D**

QR  Quick release

MS  Metal slide

D  Double D ring

MSS  1,5" metal slide buckle



## D ▶ Bed attachment

\_\_\_  Quick release ties attachment; if you choose this bed attachment, **leave blank** and go to **E**

A1  Quick release buckle 1 piece

A  Quick release 2 pieces



A1

Buckle is fixed on one strap. Offering a 10" adjustment



A

Male piece and female piece of buckle are on separate straps permitting exact adjustment.

## E ▶ Length of strap

\_\_\_  Standard , 40"; if you choose this length, **leave blank**

\_\_\_ "  choose length